



Notes: Protect the catheter and loop of tubing with an I.V. House UltraDressing® or I.V. House UltraDome®. Check IV site hourly for redness, swelling, and tenderness. **Touch, Look, and Compare** the IV site with the opposite palm and forearm. Change the TLC Wrist Splint each time a new IV is inserted.



Video Directions: www.ivhouse.com/how

Manufactured for I.V. House, Inc.
 2072 Concourse Dr.
 St. Louis, MO 63146 USA
 IVinfo@ivhouse.com (E-mail)
www.ivhouse.com
 800-530-0400 (Toll Free)
 314-453-9200 (Phone)
 314-453-9576 (Fax)

Patents Pending
 All Rights Reserved
 Made in USA



WRIST / CHILD / 12-41 KG

TLC® Wrist Splint with Straps
 Touch/Look/Compare



939L-W-Ultra

1 Large (Child) Size Foam Padded Wrist Splint with Three Straps and Velcro® Closure

Fits most children ages 3–10 years, or weighing 12–41 kg

CATALOG ITEM:

REF 939L-W-Ultra



LOT 2022-07-01-939L-W-Ultra REV 0722



I.V. HOUSE®
 Protection Over and Above®

